

Instructions for cleaning and sterilization

Refraction device
EYE REFRACT®

VISIONIX

This instruction guide aims to clean and sterilize your EYE REFRACT® for common use .
For other instructions, please refer to the USER MANUAL.

For your comfort it's recommended to clean the device regularly. Please clean the instrument following these instructions. Cleaning conducted by a process other than the one described in this instruction may cause deterioration or damage to the equipment.



1- Cleaning and sterilization of chin rest and forehead rest support

Cleaning of the chin rest and headrest should be done before each measurement. The patient's chin and head should be on the chin rest and forehead during all diagnostic procedures. That's why it's important to clean the chin rest and forehead support with Cidalkan wipes incl. alcohol. disinfectant or other alternative.


Optional items

Chin rest paper (REF. 4100-0014-00)

When using the chin rest paper, use a new sheet for each patient.



For prevention of cross-infection, replace with a new paper for every new patient. Pull and replace with a new sheet of paper.



2- Cleaning of Removable Chick protection

Chick protection are removable spare parts, it is recommended to clean the front surface between each patient and to remove them totally between each 5 patients to clean them completely.



During the front surface cleaning of Chick protection take care of the apertures and the lenses inside, do not clean do deep you may damage the lenses.

3- Cleaning of tablet and near vision screen

The tablet should not be seen or touched by the patient. But clean the screen regularly if possible between each patient, after you take off your gloves. Use Cidalkan wipes or other alternative.



Put the tablet on a flat surface previously cleaned. The tablet must be wiped gently without applying excessive force.

4- Cleaning of the internal body of the device

Between the left and right eyes, proceed to a regular cleaning of the interior of the device.

The patient will touch the device with his/her nose during the exam. Before cleaning the unit, turn it off and unplug it from the electric outlet.

To clean the plastic surface of the device gently wipe the top, bottom, and front surfaces and insist on the contact surfaces.



Only use the below recommended solutions. Other solutions may damage or discolour the device.



Avoid to give the printed results to your patients.



Wipes - Instructions for use please read carefully !

01

Take a wipe out of its packaging and close the packaging tightly after use.

02



Clean thoroughly by hand wiping. Leave the product on for 5 minutes. Do not rinse.

03

Throw away the wipe. Change wipe for each medical surface or device.

REF. 369113 Cidalkan wipes incl. alcohol disinfectant

Standards CE 0459



Luneau Technology

LUNEAU TECHNOLOGY OPERATIONS

2 Rue Roger Bonnet-27340 PONT DE L'ARCHE - FRANCE - Tel. + 33 232 989 132 • contact@luneautech.com • www.luneautech.com