

# Eye Refract™

## AI powered physiological refraction

Visionix® Eye Refract™, manufactured by Luneau Technology group, is the only AI powered device that offers « physiological refraction » to establish the most comfortable prescription for your patients.



### PHYSIOLOGICAL REFRACTION

We speak of physiological refraction because the measure is taken dynamically until a stable measure is reached: it is an adaptive process taking into account the physiological reaction of the patient to the proposed prescription binocularly, in the same way people are naturally using their eyes.

Visionix pioneered the usage of wavefront in eyecare: Eye Refract™ uses aberrometry for an accurate measure.



### ARTIFICIAL INTELLIGENCE

In addition to the physiological refraction process which is AI powered, Quick Pro algorithm in Eye Refract simulates the steps an experienced optometrist would take to converge quickly to a comfortable prescription

- It was continuously optimized analyzing hundreds of cases of subjective refractions
- Comparing results from ER versus classical subjective
- involving top Professors in Optometry & multiple experienced Optometrists

Contact us : [www.luneautech.com](http://www.luneautech.com)